

Add an Anaesthetic Alert in TrakCare

Calum McDonald
Consultant Anaesthetist
March 2021

The logo for Aberdeen Anaesthesia, featuring a stylized blue 'A' shape that frames the text 'Aberdeen Anaesthesia'.

www.aberdeenanaesthesia.org



@Abdn_Anaesth



Start

- The guide will show you how to add an anaesthetic alert e.g. airway management problems.
- Adding an allergy is a similar process, but is not covered in this presentation.
- The example in this presentation is just an example and does not mean that all grade 2 laryngoscopies should be recorded.
- The decision as to what is an appropriate alert to record is down to the individual anaesthetist's clinical judgement.

Open the EPR for your patient

The screenshot shows the TrakCare EPR system interface. The top navigation bar includes 'Home | Tools | Messages | Logout' and 'User: mcdonc2 | Location: Aberdeen Royal Infirmary | Lock'. The main header displays 'EPR' and 'Patient Find > Patient List > EPR Chartbook'. The left sidebar lists various medical categories: Allergies / alerts, All investigations, Grouped lab results, Laboratory, Radiology, Other investigations, Service requests, Observations, Clinical assessment, Documents (sdi store), Primary care summary, Historic / sdi store results, C Cube, Admin letters, Diagnosis / procedures, and Patient tree. The main content area shows the 'Allergies' and 'Alerts' sections. The 'Allergies' section has a 'New' button and a table with columns: Category, Allergen, Free Text Allergy, Nature Of Reaction, Severity, Onset Date, Comments, Update User, Update Date, Status, and Edit. The 'Alerts' section also has a 'New' button and a table with columns: Category, Alert, Message, Date Entered, Review Date, Exp. Date, Update User, Update Date, Status, Edit, and Created User. A blue arrow points from the 'New' button in the Alerts section to the right-hand text.

- Allergies / Alerts is the first page to display.
- Allergies and alerts are shown here.
- Click on New in the Alerts section to record an anaesthetic alert or hazard.
- A new window will open.

Alert Window

SCGR-T2014-LIVE - Internet Explorer

Home | Tools | Messages | Logout

User: mcdonc2 | Location: Aberdeen Royal Infirmary | Lock | TrakCare

CHI: [REDACTED] Name: [REDACTED] Gender: Male DOB: [REDACTED] Age: [REDACTED] Yrs Phone: [REDACTED]

Permanent Alerts

Please remember any User can view permanent alerts

Alert Category

Alert

Date Created 20/06/2017

Extra Information

Closed Flag ☐

Display Alert ☐

Expected Review Date

Expiry Date

User Name mcdonc2

Password

Update

- Click on the spyglass to select Alert Category and choose **Clinical**.
- Click on the Alert spyglass and choose **Anaesthetic Risk**.
- The date will default to today.

Alert Window

SCGR-T2014-LIVE - Internet Explorer

Home | Tools | Messages | Logout

User: mcdonc2 | Location: Aberdeen Royal Infirmary | Lock | TrakCare

CH: [REDACTED] Name: [REDACTED] Gender: Male DOB: [REDACTED] Age: [REDACTED] Yrs Phone: [REDACTED]

Permanent Alerts

Please remember any User can view permanent alerts

Alert Category: Clinical

Alert: Anaesthetic Risk

Date Created: 20/06/2017

Extra Information: Grade 2 laryngoscopy. Bougie and BURP used.

Closed Flag ☐

Display Alert ☐

Expected Review Date: [REDACTED]

Expiry Date: [REDACTED]

User Name: mcdonc2

Password: [REDACTED]

Update

- Enter brief details of the risk here .
- Enter your password here.
- Click Update.

Alert Added

SCGR-T2014-LIVE - Internet Explorer

Home | Tools | Messages | Logout

User: mcdonc2 | Location: Aberdeen Royal Infirmary | Lock

TrakCare by InterSystems

EPR

Patient Find > Patient List > EPR Chartbook

New Request Episode List Personal Patient List Correspondence-Referral Wards Inpatients

My Recent Patients Outpatient Worklist Schedule Order Comms Applications BOXI Reports Investigations

CHI: [REDACTED] Name: [REDACTED] Gender: [REDACTED] DOB: [REDACTED] Age: [REDACTED] Yrs Phone: [REDACTED]

Allergies / alerts

Allergies + New

Category	Allergen	Free Text Allergy	Nature Of Reaction	Severity	Onset Date	Comments	Update User	Update Date	Status	Edit
----------	----------	-------------------	--------------------	----------	------------	----------	-------------	-------------	--------	------

Alerts + New

Category	Alert	Message	Date Entered	Review Date	Expiry Date	Update User	Update Date	Status	Edit	Created User
Clinical	Anaesthetic Risk	Grade 2 laryngoscopy. Bougie a	20/06/2017			Calum McDonald	20/06/2017	Active		Calum McDonald

Documents (sc store)

Primary care summary

Historic / sc store results

C Cube

Admin letters

Diagnosis / procedures

Patient tree

mcdonc2 200%

- The alert now shows on this screen.
- There is also an Alert icon in the patient banner.
- Click on this arrow to show the full message.

The end

- The procedure for adding allergies is similar. Just click on the New link in the Allergies section.
- Any problems, contact the IT Service Desk online or by phone – 54444, or email me:
 - calum.mcdonald2@nhs.scot